

# Q1 - Model fine tuning

Fine-tuning a YOLO model for detecting medical anomalies in X-ray images means adapting an already trained YOLO model to perform well on a new task—identifying abnormalities in X-ray images. Here's a detailed and beginner-friendly explanation:

---

## What is YOLO?

- YOLO (You Only Look Once) is a popular machine learning model used for object detection.
- It takes an image as input and identifies objects in the image by drawing boxes around them and assigning labels to each box.
- Imagine you're identifying apples, bananas, and oranges in a photo. YOLO finds each fruit, draws a box around it, and says what it is.

When we use YOLO for medical images, the goal is to identify medical issues (like fractures or tumors) in X-rays instead of apples and bananas.

---

## Step 1: Understand the Pretrained Model

- A pretrained YOLO model has already learned to detect objects in general images (like animals, cars, or people) using a dataset called COCO.
  - But X-ray images are different. Medical anomalies are subtle and not present in COCO data. That's why we need to *fine-tune* the model so it can detect these specific abnormalities.
- 

## Step 2: Gather and Prepare Your Data

1. **Collect Images:** Get lots of X-ray images. Some might show normal X-rays, while others might show abnormalities.
  2. **Annotate the Images:**
    - Draw boxes around the abnormalities in the images.
    - Label each box with the type of abnormality (e.g., "fracture" or "tumor").
    - There are tools like **Labellmg** to help you draw these boxes on your computer.
  3. **Organize the Data:**
    - Save the annotations in a format YOLO understands (e.g., `.txt` files alongside each image).
  4. **Resize the Images:**
    - YOLO needs images of a specific size (e.g., 416x416 pixels). This resizing makes sure all images are the same size for the model.
  5. **Augment the Data:**
    - Data augmentation means creating variations of your images to teach the model better. For example:
      - Flip the X-ray horizontally or vertically.
      - Slightly rotate the image.
      - Adjust brightness and contrast to simulate different scanning conditions.
- 

## Step 3: Modify the YOLO Model

- The pretrained YOLO model has a part that outputs labels like "cat" or "car."
  - Replace this part with one that outputs your labels (like "fracture" or "tumor").
  - This is like swapping the signboard of a store: Instead of selling fruits, the store now sells medical supplies.
- 

## Step 4: Train the Model

1. **What is Training?**
  - Training is when the model "looks" at your annotated X-rays, learns what the anomalies look like, and adjusts itself to make better predictions.
2. **How to Start?**
  - Use your prepared dataset.
  - Load the pretrained YOLO model, but tell it to focus on learning only about your data.
3. **Settings to Adjust (Hyperparameters):**
  - **Learning Rate:** Controls how fast the model updates itself while learning. Start with a small value like 0.0001 to avoid sudden wrong updates.
  - **Batch Size:** This is how many images the model looks at at once during training. A batch size of 16 or 32 is a good starting point.
  - **Number of Epochs:** One epoch means the model has looked at all the images in your dataset once. You might need 50–100 epochs for this task.
4. **Use Validation Data:**
  - Keep some X-rays aside (don't use them for training). After training, test the model on these images to see how well it performs.
5. **Monitor the Model's Progress:**
  - During training, the model calculates "loss," which tells you how wrong its predictions are.
  - A decreasing loss means the model is learning.

---

## Step 5: Evaluate and Improve

1. **Test the Model:**
    - Use new X-rays to test the model's performance. Check:
      - Does it detect all abnormalities?
      - Does it avoid marking normal areas as abnormal?
  2. **Metrics to Measure Performance:**
    - **Precision:** How many of the detected abnormalities are correct?
    - **Recall:** How many of the actual abnormalities were detected?
    - **mAP (mean Average Precision):** A summary score that tells you how well the model is doing overall.
  3. **Tweak Settings if Needed:**
    - Adjust the learning rate or batch size.
    - Train for more epochs if it's not learning enough or stop early if it's overfitting (memorizing instead of generalizing).
- 

## Step 6: Deploy the Model

1. **Set a Confidence Threshold:**
    - The model predicts anomalies with a "confidence score" (e.g., 85% sure it's a tumor). Decide a threshold (e.g., 50%) to filter predictions.
  2. **Run the Model on Real Data:**
    - Feed new X-rays into the model. It will highlight areas with abnormalities and label them.
  3. **Validate with Doctors:**
    - Have medical professionals check the results to ensure accuracy.
- 

## Important Tips

- **Patience:** Fine-tuning can take hours or even days, depending on your hardware.
  - **Hardware:** Use a powerful computer with a GPU (Graphics Processing Unit). GPUs make training much faster.
  - **Learning from Failures:** If the model misses abnormalities or makes many mistakes, look at the images it got wrong. This can give clues about what's missing in the training data.
- 

By following these steps, you're teaching the YOLO model to "look" at X-rays and identify medical problems as a trained doctor would. Fine-tuning is a process of adapting a general-purpose tool to a specific and important task, like diagnosing medical conditions.